

CXCR4-chemokine (C-X-C motif) receptor 4 Gene

Reference: PTXBC020968

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD184; D2S201E; FB22; HM89; HSY3RR; LAP-3; LAP3; LCR1; LESTR; NPY3R; NPYR; NPYRL; NPYY3R; WHIM; WHIMS; C-X-C chemokine receptor type 4; CD184 antigen; LPS-associated protein 3; SDF-1 receptor; chemokine (C-X-C motif) receptor 4; chemokine receptor; fusin;		

For research use only.

GENE DESCRIPTION	chemokine (C-X-C motif) receptor 4
GENE ID	7852
ORIGIN SPECIES	Human
ACCESSIONS	BC020968
VECTOR	pDONR223
NCBI SYMBOL	CXCR4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggggacagatatacacttcagataacacacagaggaatgggctcaggggactatgactccatgaaggaaacctgttccgtgaagaaatgctaattcaataaaattcttccgccaccatctactccatctcttaactggcattgtgggcaatggattggatcctgtgca
tgggttacagaagaactgagaagcatgacggacaglacaggctgcacctgtcagtgccgacgtcctcttctgcatcacgcttccctctctgggagttgatgccgtgggcaaaactggatcttgggaacttctatgcaaggcagtcctatgtcatctacacagtcacacctctacagcagttg
ccatcctctggccttcacagtcctgacggcctacgtggccatgtccagccaccaacagtcagaggccaaggaagctgttggctgaaaagtggtctatgttggcgtctggatccctgcccctctgctgactattcccgaacttcttggcaacgtcagtgaggcagatgacagatatct
gtgaccgcttctaccccaatgactgtgggtgtgtgtccagtttcagcacatggttggccttactcctgctgattgtcatcctgtcctgctattgcattatcatctccaagctgtcacactccaaggccaccagaagcgaaggccctcaagaccacagtcacctcctatcctgttct
tcgctgttggctgccttactacattggatcagcagtcgactcctcctcctctggaaatcatcaagcaagggtgtgagtttgagaacactgtgcacaagtggaattccatcaccgagggccctagcttcttccactgtgtctgaacccatcctctatgcttcttgagccaaatttaaacctc
tgccagcagcactcactcctctgtgagcagagggtccagcctcaagatcctctccaaaggaaagcgaggtggacattcatctgttccactgagctgtgagcttcaagtttcaactcagctaa