



Reference: PTXBC021098

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BCL11A B-cell CLL/lymphoma 11A (zinc finger protein) isoform 1; BCL11a-M; BCL11A-XL; BCL11A-S; BCL11A-L; CTIP1; DILOS; EVI9; HBFQTL5; ZNF856; B-cell lymphoma/leukemia 11A; B-cell CLL/lymphoma 11A (zinc finger protein) isoform 2; BCL-11A; C2H2-type zinc fi		

For research use only.

GENE DESCRIPTION	B-cell CLL/lymphoma 11A (zinc finger protein)
GENE ID	53335
ORIGIN SPECIES	Human
ACCESSIONS	BC021098
VECTOR	pDONR223
NCBI SYMBOL	BCL11A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr