

APLN-apelin Gene

Reference: PTXBC021104

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	APEL; XNPEP2; AGTRL1 ligand; APJ endogenous ligand		

For research use only.

GENE DESCRIPTION	apelin
GENE ID	8862
ORIGIN SPECIES	Human
ACCESSIONS	BC021104
VECTOR	pDONR223
NCBI SYMBOL	APLN
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccgccacctggcagcccagagggtcaaggaatggccaggggccctggcaggagggtcgaggaaatcgcgcgccagcggcccgccctcccataaggacccatgcttctgaagcaggactgaaggggccccaagtgccacccccggcggttatgtctctcca
tagattggctgtctctctggaggccctcacgtccattcagctctcacctgcacctgtctgagccaccagtgggccagctctctcacctgctctcccccagtgccgtgctctggctgtagttggatgattcccttctctcacaagaatccgtccagtcacatctctggccctccctgg
actga