

MRPS16-mitochondrial ribosomal protein S16 Gene

Reference: PTXBC021106

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CGI-132; COXPD2; MRP-S16; RPMS16; 28S ribosomal protein S16, mitochondrial; S16mt; mitochondrial ribosomal protein S16		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein S16
GENE ID	51021
ORIGIN SPECIES	Human
ACCESSIONS	BC021106
VECTOR	pDONR223
NCBI SYMBOL	MRPS16
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccacctcactctcctctgcaaggcctaccgtggggccacttaaccatccgcttgcctgggtggctgcaccaatcgccgttctaccgcattgtggctgctcacaacaagtgtcccagggtggcgtttcgtagcagctgggtcctatgatccattgcccaacagtcattg
gagaaaaactcgttgcctcaacctagacaggatccgtcattggattggctgcggggccacctctctaagcctatggaaaagcttctgggtcttgccttttccctctgcacctatgatgcacaaatctgagagactgcgaaggaaacgggcacgtgaagtcctgttagcttctcag
aaaacagatgcagaagctacagatacagaggctacagaaacataa