

C15orf48-chromosome 15 open reading frame 48 Gene

Reference: PTXBC021173

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FOAP-11; NMES1; normal mucosa of esophagus-specific gene 1 protein; normal mucosa of esophagus specific 1; chromosome 15 open reading frame 48		

For research use only.

GENE DESCRIPTION	chromosome 15 open reading frame 48
GENE ID	84419
ORIGIN SPECIES	Human
ACCESSIONS	BC021173
VECTOR	pDONR223
NCBI SYMBOL	C15orf48
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagcttttccaactcctgatgaaaaggaaggactcattcccttggtggtgttcagctgtggcggcggtggagcctcatcttgcgtgtgtattctttggaaaaccgatgtgatccttgatcgaaaaaaatccagaaccttgggaaactgtggacacctgtacctcaaaagcttat
aacaatcaaccaacaatggaaacccattgaagagttgcaaatgtccaaagggtgaccaaata