

HMGB4-high-mobility group box 4 Gene

Reference: PTXBC021180

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	dJ1007G16.5; high mobility group protein B4; HMG2 like; high mobility group box 4		
For research use only.			

GENE DESCRIPTION	high-mobility group box 4
GENE ID	127540
ORIGIN SPECIES	Human
ACCESSIONS	BC021180
VECTOR	pDONR223
NCBI SYMBOL	HMGB4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggaaaaaataccagctaaagcctaagcgaatgtctcttctacgttcactttttgctgaattacagaacaaattcaaggagcagcagccaagtacctacgttggtcttaagagttctctagaaagtgttcgaaaaatggagatccatctcaagcatgaaaaggccaaatatgaagc
cctggccaaactcgacaaagcccataccagggaagaatgatgaattatgttggaagaggagaacggagaagcgggatccccaggcaccagacggcctccatcatccttctactcttctgccaagaccactatgctcagctgaagagggagaacccgaactggctgggtg
cagggtggccaaggccacagggaagatgtgtcaacagcgacagacctggagaagcacccttatgagcaagagtggtctctctgagagctaaagtacttcgaggaacttgaactctacgtaaaccaatgtaatgccaggaagaagtaccgaatgtcagctagaaccgggtgcagagggga
aaagagtcagcgagagctga