

CST1-cystatin SN Gene

Reference: PTXBC021225

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	cystatin-SN; cystatin-SA-I; cystatin 1; cystatin SA-I; cysteine proteinase inhibitor, type 2 family; salivary cystatin-SA-1; cystatin SN		

For research use only.

GENE DESCRIPTION	cystatin SN
GENE ID	1469
ORIGIN SPECIES	Human
ACCESSIONS	BC021225
VECTOR	pDONR223
NCBI SYMBOL	CST1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcccagtagtaccctgtgtcctctgtgtccaccctagctgtggccctggcctggagcccaaggaggaggataggataatccgggtggcatctataacgagacctcaatgatgagtggtacagcgtgcccttcacttcgccatcagcaggtataacaaggccacaaa
gatgactactacagacgtccgtgctgctgactaagagccaggcaacagaccgttgggggtgaattacttcttcgacgtagaggtgggcccgaaccatgtaccaagtcccagcccaacttggacacctgtgccttcgaacagccagaactgcagaagaacagttgtgctcttgc
agatctacgaagtccctgggagacagaaggctccctggtgaaatccagggtgcaagaatcctag