

KRT81-keratin 81 Gene

Reference: PTXBC021241

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HB1; Hb-1; KRTHB1; MLN137; ghHkb1; hHAKB2-1; keratin, type II cuticular Hb1; K81; MLN 137; ghHb1; hair keratin K2.9; hard keratin, type II, 1; keratin 81, type II; keratin, hair, basic, 1; metastatic lymph node 137 gene protein; type II hair keratin Hb1;		

For research use only.

GENE DESCRIPTION	keratin 81
GENE ID	3887
ORIGIN SPECIES	Human
ACCESSIONS	BC021241
VECTOR	pDONR223
NCBI SYMBOL	KRT81
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagggccacgggtgatcaggcacggggagaccctgcggcgcaccaaggaggagatcaacgagctgaaccgcatgaccagaggctgacggccgaggtggagaatgccaagtgcagaaactccaagctggaggccggtggcccagctctgagcagcagggtgaggcggccc
tcagtgatccccgtgcaagctggccgagctggaggcgccctgcagaaggccaagcaggacatggcctgcctgacagggagtagaccagggagtgatgaactccaagctggccctggacatcgagatcgccactacaggcgctgctggaggcgaggagcagaggctatgtga
aggcattggggctgtgaatgtctgtgtcagcagctccccggggcggtgtgtgcggggacacctgtcgtgtcaggctccccggccagtgactggcagtgctgtcagcgcctccgtgcaacgggaacgtggcgggtgagcaccggcctgtgtgcgcctgcggccaattgaacaccacctgc
ggaggggggttctcggcgtgggctcctgtgtatcagctccctgggtgtgggtcttcgggcagcagctgccggaatgttag