



Reference: PTXBC021252

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	polycomb protein SCMH1; Scml3; sex comb on midleg homolog 1 (Drosophila)		
For research use only.			

GENE DESCRIPTION	sex comb on midleg homolog 1 (Drosophila)
GENE ID	22955
ORIGIN SPECIES	Human
ACCESSIONS	BC021252
VECTOR	pDONR223
NCBI SYMBOL	SCMH1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr