

AIF1L-allograft inflammatory factor 1-like Gene

Reference: PTXBC021253

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	2810003C17Rik; AI043124; C87647; Iba2; allograft inflammatory factor 1-like; ionized calcium binding adapter molecule 2		

For research use only.

GENE DESCRIPTION	allograft inflammatory factor 1-like
GENE ID	83543
ORIGIN SPECIES	Human
ACCESSIONS	BC021253
VECTOR	pDONR223
NCBI SYMBOL	AIF1L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcggcgagctcagcaacaggtccaaggagggaaggcgttcgctgctcaagccggcaggagaggcgtggccgagatcaaccgggagttctgtgtgaccagaagtacagtgatgaagagaacctccagaaaagctcacagcctcaagagaagtacatggagttg
acctgaacaatgaaggcgagattgacctgatgtctttaagaggatgatggagaagcttgggtgtccccaagaccacctggagatgaagaagatgatctcagaggtgacaggaggagtcattgtttga