

C9orf116-chromosome 9 open reading frame 116 Gene

Reference: PTXBC021261

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0691 protein C9orf116; PIERCE1; RbEST47; p53-induced expression 1 in Rb-/- cells; p53-induced expression in RB-null cells protein 1; pierce 1; chromosome 9 open reading frame 116		

For research use only.

GENE DESCRIPTION	chromosome 9 open reading frame 116
GENE ID	138162
ORIGIN SPECIES	Human
ACCESSIONS	BC021261
VECTOR	pDONR223
NCBI SYMBOL	C9orf116
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgaggaaatgccccagagcgtgcgcggagcctgtggcgcccaaggccacggccccccggagagagaccagcgactaccgctgagcgcggacctgccggcagggttaacaacccgggggtgtccgggctacaggaccagaaggctgtctccgtgtacaggac
cagtaaccaggccttacgggagcagagccccaccgtgcacgagatgcctcggtgggaatgtccgggaacaatactctcaatgtttacctggagaaaagcatcgtga