

MGC29506-hypothetical protein MGC29506 Gene

Reference: PTXBC021275

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical protein MGC29506
GENE ID	51237
ORIGIN SPECIES	Human
ACCESSIONS	BC021275
VECTOR	pDONR223
NCBI SYMBOL	MGC29506
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggctgtcactgccactgctgctgctgctgggagcctgggcatccagggggcctcggggacaggggcgcactcacagccacagcccaactggatgaggagatgactcagccacatgccgctcactgcgctgtgatgcctgcagagctgtggctaccaga
tftggcaaatctggcaaggcagagaccaactcatactcaaaactctggggcgggcgagctgagcagttggtctacagggatgtcctggaccggagctgctccgggaactggcaggactacggagttcgagaagtggaacagtgaaactctcacagggcaggactta
gcgagggccagagccaagcagcgtgatgtcacagggggccctggcctaccaggtctccaggacatgtttgcactactgggggagtttgagaagaccagatctatgaagccaccaacaaggccgaggggctctggagcattgctatgtggggacccagggggcc
tgctcagagaaggtgtcagccacaagagaagagctctag