



Reference: PTXBC021571

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	E3 ubiquitin-protein ligase RNF139; HRCA1; RCA1; multiple membrane spanning receptor TRC8; patched related protein translocated in renal cancer; translocation in renal carcinoma on chromosome 8 protein; ring finger protein 139		

For research use only.

GENE DESCRIPTION	ring finger protein 139
GENE ID	11236
ORIGIN SPECIES	Human
ACCESSIONS	BC021571
VECTOR	pDONR223
NCBI SYMBOL	RNF139
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr