



Reference: PTXBC021578

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CIC-2; CLC2; ECA2; ECA3; EGI11; EGI3; EGMA; EJM6; EJM8; LKPAT; cIC-2; chloride channel protein 2; chloride channel 2; chloride channel, voltage-sensitive 2; chloride voltage-gated channel 2		

For research use only.

GENE DESCRIPTION	chloride channel 2
GENE ID	1181
ORIGIN SPECIES	Human
ACCESSIONS	BC021578
VECTOR	pDONR223
NCBI SYMBOL	CLCN2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]