

NAIF1-nuclear apoptosis inducing factor 1 Gene

Reference: PTXBC021580

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C9orf90; bA379C10.2; nuclear apoptosis-inducing factor 1		
For research use only.			

GENE DESCRIPTION	nuclear apoptosis inducing factor 1
GENE ID	203245
ORIGIN SPECIES	Human
ACCESSIONS	BC021580
VECTOR	pDONR223
NCBI SYMBOL	NAIF1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccgtcccagcaagaaaggaagatgaactctcagagcgggaggtggagatcatcgtggaggagctggagctgaagaagcacctgtggtgaaccactcaacgccgggtaccctggccgccaagatgcggcctggcacggcatcctgagaagggtcaacgccgtggccacctgcccagagagctgctgaggtcaagaagaagtgtctgacctcaagaccgaggctccgtcgaagggtgccaggtccggccgctggagggtggtgaggcgcggggccacctgaggaggacggagctgggggctgggacaggcgtggcagtggtggcgggtggcccagctgtagccccagtgctgtgacccccatgaacaacgtatctgcaacctgtgtggcgaggccaccatcatcagctgccagaccacagagatccacctgtggccctcggaacctcgggcaccgcagccgcagccacgtcacctgacacagagtgtgta