

STK32A-serine/threonine kinase 32A Gene

Reference: PTXBC021666

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	YANK1; serine/threonine-protein kinase 32A; A930015B13Rik; CTB-108O6.2; yet another novel kinase 1; serine/threonine kinase 32A		

For research use only.

GENE DESCRIPTION	serine/threonine kinase 32A
GENE ID	202374
ORIGIN SPECIES	Human
ACCESSIONS	BC021666
VECTOR	pDONR223
NCBI SYMBOL	STK32A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggagccaacactcaagaaaccaccagtgtttgatgaaatgaagatgtcaactttgaccactttgaaattttcgagccattgggaaggcagttttgggaaggctgcattgtacagaagaatgataccaagaagatgtacgcaatgaagtacatgaataaacaagaagtcgtggag
cgcaatgaagtgaagaatgtcttcaaggaaactccagatcatgcagggctgtgagcaccctttctgttaatttgggtattcctccaagatgaggaagacatgttcattgtgtgggacctcctctgtgggtgagacctgcgttatcacctgcaacagaacgtccactcaaggagaacag
tgaagctctcatctgtgagctgtcatggccctggactacctgcagaaccagcgcatcattcacagggatgaagcctgacaataatttactgacgaacatgatacctggctctctacaagtccactga