

ANKRD22-ankyrin repeat domain 22 Gene

Reference: PTXBC021671

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ankyrin repeat domain-containing protein 22; ankyrin repeat domain 22		
For research use only.			

GENE DESCRIPTION	ankyrin repeat domain 22
GENE ID	118932
ORIGIN SPECIES	Human
ACCESSIONS	BC021671
VECTOR	pDONR223
NCBI SYMBOL	ANKRD22
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggaatcctatactctgagcccatctgccaagcagcctatcagaatgactttggacaagtgtggcggtgggtgaaagaagacagcagctatgccaacgttcaagatggctttaatggagacacgccctgatctgtcttgaggcgagggcatgtgagaatcgttctctctttaagaa
gaaatgctaattgtcaacctcaaaaccagaagagagaacctgctgcatcatgctgtgaagaaaaattaccctcattgattatctactaattatcctcttaattgctgtctgtcttattgggtatttctcatggtatcaagacaaagcagaatgaggctctgtacgaatgctactgtgctggt
gtcgaagttaattgtacagattgttatggctgtaccgattacattatgacctgtgaaatgaaaaccagctcttattccctctgctcttggagcccggtgagaccccaataaagaataagcatggtgagagctcactggatattgcacggagattaaaatttccagattgaattatgctaag
gaaagcattgtaa