

C5orf46-chromosome 5 open reading frame 46 Gene

Reference: PTXBC021680

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C5orf46; SSSP1; CTC-327F10.5; skin and saliva secreted protein 1; chromosome 5 open reading frame 46		

For research use only.

GENE DESCRIPTION	chromosome 5 open reading frame 46
GENE ID	389336
ORIGIN SPECIES	Human
ACCESSIONS	BC021680
VECTOR	pDONR223
NCBI SYMBOL	C5orf46
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgtcttagtacttcgcctgacagttgtcctgggactgttcttattctgacctgctatgcagacgacaaaccagacaagccagacgacaagccagacgactcgggcaaagacccaaagccagacttcccaaatcctaagcctcctgggcacagatcattgagaatgcagt
cgagttcatcctccgtccatgtccagagcacgtaa