

C9orf61-chromosome 9 open reading frame 61 Gene

Reference: PTXBC021685

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C9orf61; X123; protein FAM189A2; Friedreich ataxia region gene X123; family with sequence similarity 189 member A2		

For research use only.

GENE DESCRIPTION	chromosome 9 open reading frame 61
GENE ID	9413
ORIGIN SPECIES	Human
ACCESSIONS	BC021685
VECTOR	pDONR223
NCBI SYMBOL	C9orf61
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaccgcaggatggttgctctgatgtattccccctgccacacatctacggagctcgaatcaaaggtgtggaagtgtctctctggtatccccccgcgcctatgaagctgtggtgagccagatggaccaggagcagggatcttcattccaaatgtcagaaggatcagaagctgctgtg
atccccatggatcgggctgcacacaagtgtactcaagatggggacattcctaatacattcgcgaagaaatgcattccacctcaactccagttcaacccctgggtgctctatcagaagccggagagccctcccaccttgaggaccaggtcgaagagtgacctgtgctccatctctga
ggagagagataatggcccatgtcagccagcacatccagggcccaagctgccctcgggagacagcctggcctgctgcacctccagagctgcggcgacctcacaccttcacaccagcggggaggccccgagccgagaggagggccccggcgagtgaggctga