



Reference: PTXBC021708

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HLA-G histocompatibility antigen, class I, G; MHC-G; HLA class I histocompatibility antigen, alpha chain G; HLA G antigen; MHC class I antigen G; b2 microglobulin; major histocompatibility complex, class I, G		

For research use only.

GENE DESCRIPTION	major histocompatibility complex, class I, G
GENE ID	3135
ORIGIN SPECIES	Human
ACCESSIONS	BC021708
VECTOR	pDONR223
NCBI SYMBOL	HLA-G
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]