

SRrp35-serine-arginine repressor protein (35 kDa) Gene

Reference: PTXBC021715

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SRrp35; SFRS13B; SFRS19; serine/arginine-rich splicing factor 12; 35 kDa SR repressor protein; SR splicing factor 12; serine-arginine repressor protein (35 kDa); splicing factor, arginine/serine-rich 13B; splicing factor, arginine/serine-rich 19; serine a		

For research use only.

GENE DESCRIPTION	serine-arginine repressor protein (35 kDa)
GENE ID	135295
ORIGIN SPECIES	Human
ACCESSIONS	BC021715
VECTOR	pDONR223
NCBI SYMBOL	SRrp35
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctcgtacacgagggcccccaacacccctgttcacaggaacgtcgggacgccaccaggccctgaggacttgccggtgagtttggtcgatagggccctatagtagacgtttacattccacttgacttctacactcggccccaagaggatttgcttatgttcaatttgaagatgttcga
gggtcgtgaagatgctctttataacctcaatagaaatgggtatgtggccgtcagattgaatacagtttgcacaaggtgatcgcaaaacaccaggcccaatgaaatcaaaagaacgtcatccttcttctccaagtgatcacaggagatcaagaagccccagccaaagaagaactcgaagtag
aagttctcatggggaagaataggagcgggtcagacagcctaaagagctcgcacacagcgcgatttcttatagcaagtcctaaatcgttccaaatcattaccaagcgggtctacctcagcaaggcagtcagaactccaagaaggaaatttggctctagaggacggtcaaggtccaagtc
cttacaagaaggtccaagtcgaataggaaatcacagtcaggttcacctcaaaagcagactagctcaggaacaaatcaagatcacatggaagacattctgactcaatagcaagatccccgtgtaaatctccaagggtataccaattctgaaactaaagtacaacagcaagcatttca
tttccgtcacattccagatctcgaagttatctcataaaacagttgggtga