

FHL5-four and a half LIM domains 5 Gene

Reference: PTXBC021723

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	1700027G07Rik; ACT; FHL-5; dJ393D12.2; four and a half LIM domains protein 5; LIM protein ACT; activator of cAMP-responsive element modulator (CREM) in testis; activator of cAMP-responsive element modulator in testis; four and a half LIM domains 5		

For research use only.

GENE DESCRIPTION	four and a half LIM domains 5
GENE ID	9457
ORIGIN SPECIES	Human
ACCESSIONS	BC021723
VECTOR	pDONR223
NCBI SYMBOL	FHL5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacaactgctcattttactgtcaatactgcacagcatcattcttgggaagaataatgtactaaagatgacagtcacatactgtgttacatgttatgatcgtgtatttttaactattgcgaggaatgcaaaaaaccaattgaatctgattctaaggatctttgtacaagaccggcactggcatgaaggatgcttcaagtgcaccaaatgcaatcacttttgggaaaaagcctttgtgcgaaggatgagcgctgtgtgcacggagtgctatttcaacagagtgctctcctcaagtgcttccactgcaagaggaccatcatgctgttcccgcaaaatggaattaaaggaaactactggcatgaacctgttttgggtgtgagaattgccgacaacatagggacaagcctttgatctcaagagagtgcaattatgtgtgccaatgtttgagaaggagttgtcactactgcaacttttgaagaaggatgataactcagggtgggataacatttttgaccagcatatggcataaagagtgtttctgttagtgactgtaggaaagatctctgtgaagaacagttcatgtccagagacgactatcatttgcagtgactgctacaacatctttatgccaacagtggttagcctgttccaacccattagtggtctcacaggtgccaagtttatctgtttcaagacagcagtgccatagcgaatgctttaa