

LYZL1-lysozyme-like 1 Gene

Reference: PTXBC021730

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KAAG648; LYC2; LYZD1; PRO1278; bA534G20.1; lysozyme-like protein 1; lysozyme D1; lysozyme like 1		

For research use only.

GENE DESCRIPTION	lysozyme-like 1
GENE ID	84569
ORIGIN SPECIES	Human
ACCESSIONS	BC021730
VECTOR	pDONR223
NCBI SYMBOL	LYZL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaggacgtccctgagctgctcaccgactaggtggagcagtggtttctccgagactcaactgagaagtcagcctctggggcaggcaccaggaatctgctttcagttctgtctccggcagcgttgaggatgaaggctgaggcattctgacccctattggctgcctgtcaca
ggcgcgagtcctcaaaatctacactcgttgcaactggcaaaatattctcagggtgctggcctggacaattactgggcttcagccttggaactggaatctgcatggcatattatgagagcggctacaacaccacagccccgacggctcgtgatgacggcagcatgactatggcatctcca
gatcaacagcttcgctgggtgagacgaggaaagctgaaggagaaacactgcatgtcgcctgctcagccttgatcactgatgacctcacagatgcaattatctgtccaggaaaattgttaagagacacaaggaaatgaactattggcaaggctggaagaacattgtgaggcag
agacctgtccgagtggaagaaaggctgtgaggtttcctaa