

C10orf82-chromosome 10 open reading frame 82 Gene

Reference: PTXBC021737

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C10orf82; chromosome 10 open reading frame 82		
For research use only.			

GENE DESCRIPTION	chromosome 10 open reading frame 82
GENE ID	143379
ORIGIN SPECIES	Human
ACCESSIONS	BC021737
VECTOR	pDONR223
NCBI SYMBOL	C10orf82
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagcctccaagacctcatgagaaactgccaatcacaccaggctatagcggctttgtgccattctcagctgccaagggaatgtccaaggaggatgacatgaaccactgtgtgaaaaccttcaggagaaaacacagcgctataaagaacagctgcgggaattgtctgcgcagtgg
ccactgccccgaaactgaaacctgtcaactccgaggagacggctctgcaggcctgcaccagtacaatctgcagtaccacccctgacctggaatgcaaatatgtaagaaacctctccaggagccccgatccctggctggcaggctacctgccgagagccaaggctactgaattg
gctgtggcacgagatacactgtcatggccaaaaactgctacaaggacttctggagatcacggagaggccaagaaggcacatctgaaacatatgaagagtga