



Reference: PTXBC021825

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AVP; IFN-alpha-REC; IFNAR; IFNBR; IFRC; interferon alpha/beta receptor 1; CRF2-1; IFN-R-1; IFN-alpha/beta receptor 1; alpha-type antiviral protein; beta-type antiviral protein; cytokine receptor class-II member 1; cytokine receptor family 2 member 1; inte		

For research use only.

GENE DESCRIPTION	interferon (alpha, beta and omega) receptor 1
GENE ID	3454
ORIGIN SPECIES	Human
ACCESSIONS	BC021825
VECTOR	pDONR223
NCBI SYMBOL	IFNAR1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr