

LOC554202-hypothetical LOC554202 Gene

Reference: PTXBC021861

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical LOC554202
GENE ID	554202
ORIGIN SPECIES	Human
ACCESSIONS	BC021861
VECTOR	pDONR223
NCBI SYMBOL	LOC554202
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggattcctggaaatacctcctcaaggccagtgtagagccctcaatcaccactgccccagcttctctctccatctcctgcaccaagcaggctcaggtgtccagctgctgatgacgtaaagtgtggagttgtccaggaaggtcatgttccatgcagcagaggagcgctttgtgtg
agaagttgaagacctgctgtaacacatcatag