

LOC652968-hypothetical protein LOC652968 Gene

Reference: PTXBC021927

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical protein LOC652968
GENE ID	652968
ORIGIN SPECIES	Human
ACCESSIONS	BC021927
VECTOR	pDONR223
NCBI SYMBOL	LOC652968
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagctgcacatctagaacacccgggtgcgggtgctgagcgtcgcgcgtccgtctctgctctacacccacccgctcatcaagctgctctctgccccgcgcagccggtgcaagttctcagcctgacggagaccctgaggattacacgcttatggggacgaggagggttta
aagagctgcccccatctgagttcctgcaagtagctgaggccacatggctgctgaacgtgctgctcaccagcgggtgcggcagtgagggctgctgggtcaccagatgcccgttgcgtcatcgccactggccgagcaccaggtgctgctgatgctgtccactaccagacgg
acttcatctgtgctgggagcagacgtgctgctgctgacacacgtgcccaggagttcgaattaccgaggtggcgagagcctgctgctgacgagggatgattccagcaatggcttccccgactcagcatggggccagccccacgggtgcatcccatccagagccca
cagaaccgtctgtgctcctacactggaccctgagacgttccagccatcgccaccacccatagatgtccttctactcgcacagcaccaccaaggagcagcctctagcagtcctgaaccagctccatcagttcttgccttctccctcaggggtatctcctcattgtcatggat
gctgaacacagaaaaagttcccagtgacctcctgctgaccagctcctcggggagctgtggaggtggtgctgcgcacgtggtgacagccctgggttgatgaatgtggcatgtggcacagattgaggtccctggctgcccgtgacatctctgctactacatcagaccttcaactt
cgaccacgctgtgctggcagagcgtatcgcagcgtcaggtcctccagcggcgaggaaggcctgcttctga