

TAF11-TAF11 RNA polymerase II, TATA box binding protein (TBP)-associated factor, 28kDa Gene

Reference: PTXBC021972

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TAF11 RNA polymerase II, TATA box binding protein (TBP)-associated factor, 28kDa; MGC:15243; PRO2134; TAF2I; TAFII28; transcription initiation factor TFIID subunit 11; TAF(II)28; TAFII-28; TATA box binding protein (TBP)-associated factor, RNA polymerase I		

For research use only.

GENE DESCRIPTION	TAF11 RNA polymerase II, TATA box binding protein (TBP)-associated factor, 28kDa
GENE ID	6882
ORIGIN SPECIES	Human
ACCESSIONS	BC021972
VECTOR	pDONR223
NCBI SYMBOL	TAF11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacgatgccacgagtcgcccctccgacaaaggtggagagacaggggagtcggatgagacggccgctgtgccgggggacccgggggctaccgacaccgatggaatccagaggaaactgacggagacgcagatgtggactgaaagaagctgcagcgggaggaagcgcgagc
tcgagagtcaggatgtctcagatttaacaacagttgaaaggggaagactcattcatttaacctgcagccaaaaactgaaaatagataccaaagaaaaagagaaaaagcagaagtagatgaagatgagatcagaagatgcaaatcctgtttctcttttctgaggagcagctgaa
ccgttatgaaatgtatcgcgcgtcagctttccctaaggcagccatcaaaagcgtgacccagtcacatcggcacctctgtgtctcagaatgtttattgtctatgtctgtattccaagggttttcgtcggggaggtggtagaagaagcactggaatgtgtgtgagaagtgaggagaatgccacc
actacaacccaacatatgagggaagccgttagaaggttaagtcacaaaggacagatccctaactcgaagcacaaaaaatcatcttctctag