

NDFIP2-Nedd4 family interacting protein 2 Gene

Reference: PTXBC021988

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	N4WBP5A; NEDD4 family-interacting protein 2; MAPK-activating protein PM04 PM05 PM06 PM07; NEDD4 WW domain-binding protein 5A; NF-kappa-B-activating protein 413; Nedd4 family interacting protein 2		

For research use only.

GENE DESCRIPTION	Nedd4 family interacting protein 2
GENE ID	54602
ORIGIN SPECIES	Human
ACCESSIONS	BC021988
VECTOR	pDONR223
NCBI SYMBOL	NDFIP2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatcaccaccagccgggactgggctaccaggtgcttcttaagaggataactcagaatcatcgctatagacagccactacttcaaacccagcaccgcagattgtgcaggctgcgtcttcagcaccagcactgaaactgactctccctccaccatagtagtattact
gtggaagtacctacaactcagatcacagaagtttacgggtgagttttatcccggtgccacctccctatagcgttgctacctctcttcacatacgaatgaagctgagaaggctaaagctgctgcaatggcagctgcagcagcagaacatctcaagaattcaggaggaaggtgccaccaaga
gatgacttcagtgatgcagaccagctcagagtgagggaatgatggcattttcatgctgcatttttcatggcatttttcaactggccttggattttgtatcctctgtatcaccaataccatagctggaaggtatggtgctatctgcggtttggcctttccttgatcaaatggatcctattgtcaggtt
ttctgattttttactggatatttcaatggacagattggctttgggtgatatttctgtacttggcctgctcctttcttcagaggatttgaattatctaaaagtcagaacatgctgaaagtatggcagctgctcatagaacaaggatttcttctattgttag