

BNIP3-BCL2/adenovirus E1B 19kDa interacting protein 3 Gene

Reference: PTXBC021989

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NIP3; BCL2/adenovirus E1B 19 kDa protein-interacting protein 3; BCL2/adenovirus E1B 19kDa interacting protein 3; BCL2 interacting protein 3		

For research use only.

GENE DESCRIPTION	BCL2/adenovirus E1B 19kDa interacting protein 3
GENE ID	664
ORIGIN SPECIES	Human
ACCESSIONS	BC021989
VECTOR	pDONR223
NCBI SYMBOL	BNIP3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgcagaacggagcgcgggatgcaggagagagcctgcagggtcctggtagaactgcactcagcaataatgggaacggggcagcgtccagcctcgtttctatttataatggagacatggaaaaatactctgggacgcacagcatgagctcggacggagtagctcca
agagctctcactgtgacagccacacctcgtcgcagacaccacaagataccaacagagctctgaaacagataccatagcattggagagaaaaacagctcacagctcgggaagatgataattgaaagaaggaaagaagtgaagcatcttgaagaaaaactcagattggatgggatt
ggfcaagtcggcggaaaatattcccccaaggagttctctttaaacacccgaagcgcacggccacctcagcatgaggaaacacgagcgtcatgaagaagggggcataattctcgcagaattctgaaagtcttccatctcgtcgtctctcatttctggccatcggattggggatc
tatattggaagcgtctgacaacctcaccagcacctttga