

**PRKRIR-protein-kinase, interferon-inducible double stranded RNA dependent inhibitor,
repressor of (P58 repressor) Gene**

Reference: PTXBC021992

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PRKRIR; DAP4; P52rIPK; THAP0; 52 kDa repressor of the inhibitor of the protein kinase; 58 kDa interferon-induced protein kinase-interacting protein; P58IPK-regulatory protein; THAP domain-containing protein 0; THAP domain-containing protein 12; death-asso		

For research use only.

GENE DESCRIPTION	protein-kinase, interferon-inducible double stranded RNA dependent inhibitor, repressor of (P58 repressor)
GENE ID	5612
ORIGIN SPECIES	Human
ACCESSIONS	BC021992
VECTOR	pDONR223
NCBI SYMBOL	PRKRIR
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggacaactcaaatcaatcgtcggagggaacaccatgctgacatgtatagaagtgaactacccaatcctgacacgctgtcagctgagcttcattgttgagaatcaaatggaacacaggggaaagatatagactccgtccaccatctatgaagccctccacctgacctgacatcaag
tttttccctaatgtgatgcatgctgaaggtcctgtgtattcttctgtgatgaaggttgagaatgagcggatgaaaatggacgaaagcgtcttaagcatattgaggaaacttggacagaccaaagggtcaagtaacttggcttgccttaacataaatttgatataaacacagacctggattaat
gggtggacacatatataaactctatacaagtaagtcagagcttctacagataatccgaaactgtggaaaatacctaa