

## CRISP2-cysteine-rich secretory protein 2 Gene

Reference: PTXBC022011

|                |                                                                                                                                                                                                                                                                 |                     |                |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                                                             | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                                                            | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | CRISP-2; CT36; GAPDL5; TPX1; TSP1; cysteine-rich secretory protein 2; cancer/testis antigen 36; glyceraldehyde-3-phosphate dehydrogenase-like 5; testicular tissue protein Li 43; testis specific protein 1 (probe H4-1 p3-1); testis-specific protein TPX-1; c |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | cysteine-rich secretory protein 2                                                                 |
| <b>GENE ID</b>          | 7180                                                                                              |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC022011                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | CRISP2                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atggcttactaccggtgtgtttctgttactgtctgtccatctttacctgcagaaggaaagatcccgttttactgtttgttaaccaccagttgcaagtgcaggaggagattgtaataaacacataactaaggaaagcagctctccacctgccagtaacatgctaagatggaatg  
gagcagagaggtaacaacgaatgcccaagggtggcaacaagtgcactttacaacatagtgatccagaggaccgcaaaaccagtacaagatgtgggagaaftctctatatgtcaagtgcacctacttctgtctctgcaatccaaagctggtatgacgagatcctagattttgtctatggt  
gtaggaccaaaagagtcctaatgcagttgttgacattactcagctgtttgttactcactaccaggttaggctgtggaattgcctactgtcccaatcaagatagctctaaataactactatgtttgccaatattgtcctgctggtaataatgaaatagaagaatacccggtaccaacaaggaaac  
acctgtgccgggtgccctgatgactgtgacaaggactatgcacaaatagttgccgatcaagatctctaagtaactgtgattccttgaagaatacagctggctgtgaacatgagttactcaaggaaaagtcaaggctacttgccctatgtgagaacaaatttactga