

VSNL1-visinin-like 1 Gene

Reference: PTXBC022012

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HLP3; HPCAL3; HUVISL1; VILIP-1; visinin-like protein 1; VLP-1; hippocalcin-like protein 3; visinin like 1		

For research use only.

GENE DESCRIPTION	visinin-like 1
GENE ID	7447
ORIGIN SPECIES	Human
ACCESSIONS	BC022012
VECTOR	pDONR223
NCBI SYMBOL	VSNL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggaagcagaatagcaaaactgcccctgaagtgatggagacctggtgaagagcacagagtttaagagcatgaactcaagcagtggtacaaaggatttcaaggactgtccaagtgaggagctaaatctcgaggaattcagcagctctatgtgaattcttcttatggagacgc
ctccaagttgcccagcatgcttcgaaccttcgacaagaatgggacggcaccattgacttcgagagttcatctgcgctctgtccatcacctccagggcagctttgagcagaagctgaactgggcttcaatgtatgacctggtggtgagcgaatcacccagtgagatgct
ggagatcatcgaggctatctacaaaatgtaggcactgtgatcatgatgaaaatgaatgaggatggcctgacgcctgagcagcgagtagacaagattttcagcaagatggataagaacaagaatgaccagattacactggatgaattcaagaagctgcaaaagagcgaccttccattgtat
tactctgcagtcgacatccagaatga