



Reference: PTXBC022031

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C6orf199; AK 9; AKD1; AKD2; C6orf224; dJ70A9.1; adenylate kinase 9; adenylate kinase domain containing 1; adenylate kinase domain containing 2		

For research use only.

GENE DESCRIPTION	chromosome 6 open reading frame 199
GENE ID	221264
ORIGIN SPECIES	Human
ACCESSIONS	BC022031
VECTOR	pDONR223
NCBI SYMBOL	C6orf199
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]