

C7orf36-chromosome 7 open reading frame 36 Gene

Reference: PTXBC022043

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C7orf36; GK003; yae1 domain-containing protein 1; Yae1 domain containing 1		
For research use only.			

GENE DESCRIPTION	chromosome 7 open reading frame 36
GENE ID	57002
ORIGIN SPECIES	Human
ACCESSIONS	BC022043
VECTOR	pDONR223
NCBI SYMBOL	C7orf36
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgtgggttcaagcagcctctgatccagggccctggagacaaaggggacgtgtttaggaagaagcagacgagtcgctcctggcgagcgggaatggcagagtaaacatgcaaaagacgagtcgaagaaggttatagatggaatagatgctggcaagcagttactctcaac
agggttcaatcaaggttataagaaggtgcagaagtcatttaaaactatggagactccgaggaacattgagtgcttctcctctggttcaccttcataataaattcaactttgatcaataaaataaacaattctctggatgcagttggccagtgtgaagagatgtgctcaaacatctgaat
caatcactccacgtccatgtttagatttattggactccattgaggatagggacctttgcatgttagttcagctgagaaaaagattgatgaagctaaagatgaaagactctgtgaaataatgctggtttaacaaaactgtagcaagagccatagtgaggatgttgcataatgtagaatgtt
gtagaacacaggagcatgcacattcagaaaacccaagccccacatggattttggaacagacagccagtttagttaaacagctgggcctatcagtagatgtattacaacacctcaaacataaa