

MORN4-MORN repeat containing 4 Gene

Reference: PTXBC022054

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C10orf83; bA548K23.4; rtp; MORN repeat-containing protein 4; 44050 protein; protein 44050; retinophilin; retinophilin homolog; MORN repeat containing 4		

For research use only.

GENE DESCRIPTION	MORN repeat containing 4
GENE ID	118812
ORIGIN SPECIES	Human
ACCESSIONS	BC022054
VECTOR	pDONR223
NCBI SYMBOL	MORN4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaccctgacaaaaggttccttcacactaccagtgagggaataatcgtggcgagtggaaggagggtgagaaggaccctggggagatccatgatgaacacttcatttgctggaggacaaatacatcaggatatatag