

C8orf68-chromosome 8 open reading frame 68 Gene

Reference: PTXBC022082

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C8orf68; ERICH1 antisense RNA 1		

For research use only.

GENE DESCRIPTION	chromosome 8 open reading frame 68
GENE ID	619343
ORIGIN SPECIES	Human
ACCESSIONS	BC022082
VECTOR	pDONR223
NCBI SYMBOL	C8orf68
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggactgctcaggttttgcctggcattctgcagaagcattgctgtatctaccagacaggaatacagccgtctcctcaactatagtcctattggactgggtgcacaaatatggattgggggtgctgaattcagtcctatagcaccaccaaaccacagaagctgccccgagctgatgtcc
cgctatcggtgctggctgcaggacactgtccattgctgcatggagcctctcgtgtctgagcctcactga