

SFXN2-sideroflexin 2 Gene

Reference: PTXBC022091

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	sideroflexin-2; sideroflexin 2		
For research use only.			

GENE DESCRIPTION	sideroflexin 2
GENE ID	118980
ORIGIN SPECIES	Human
ACCESSIONS	BC022091
VECTOR	pDONR223
NCBI SYMBOL	SFXN2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggtgacgtgtctgctttaaactcgtatgcccccggtggaccagcgcaccttctggggagagtgaaagcacttctaaacatcacggacccccgactgtctttgtatctgagcgggagctggactgggccaagtgatggtggagaagagcaggatgggggtgtgccccca
ggcaccacaagtggagcagctgctgtatgccaagaagctgtatgactggccttccaccccgacactggggagaagatgaatgtatcggcgcatgtctttcagcttctggcgcatgatcacgggcttcatgctccagtctacaggacgatccggcggtgatcttctgacgtg
ggtgaaccagctcttcaatgcttgaactacaccaacagggaatgctggcttccccacatcagtcaggcagatggcccttcttcttccacagccacaaccactgctgtggccacggctgtggcgatgaacatgttgacaaagaaagcgcgcccttgggtggccgctgggtgcccccttg
ccgctgtggctgctggcctaactgtgtcaatatccccatgacgacagcaggagctcataaagggaatctgctgtaaggacaggaatgaaatgagattgttattcccgagagctgcccataggcatcacccaagtagttatttctcgatcacatgacgtcctgggatgatcttgc
tgccagtcacatggaagccttgagaaattgcacttcacagaaagtaaggtcctgcacgccccattgcaggctatgctgagcgggtgcttctcatcttcatggtccagtggcggtgtgggttttccacagaaatgtaattgccagtttctctatctggaaccgaagctccaagacac
atcaaggccaagtatggagaacttgagccttatgtctacttcaataagggtctctaa