



Reference: PTXBC022243

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	methylcytosine dioxygenase TET3; hCG_40738; tet oncogene family member 3; tet methylcytosine dioxygenase 3		

For research use only.

GENE DESCRIPTION	tet oncogene family member 3
GENE ID	200424
ORIGIN SPECIES	Human
ACCESSIONS	BC022243
VECTOR	pDONR223
NCBI SYMBOL	TET3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]