



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

TRA@-T cell receptor alpha locus Gene

Reference: PTXBC022317

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

TRA@; IMD7; TCRA; TCRD; TRAC; T-cell antigen receptor, alpha polypeptide; T-cell receptor, alpha (V,D,J,C); T-cell receptor alpha locus

For research use only.

GENE DESCRIPTION	T cell receptor alpha locus
GENE ID	6955
ORIGIN SPECIES	Human
ACCESSIONS	BC022317
VECTOR	pDONR223
NCBI SYMBOL	TRA@
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgtttcttcacgctctgtgtgtatttttggcccttcagctcatctgcatcaagctgtggcccaagaaggttactcaaacccagctcatcagtatcatgcagctgaggagaagacgctaccctgaactcgtctgatgaacaagttgtgtgtcatatfatatatttttggtaacaagaacttcccgcaa
agagatgattttccttattcggcagggttctgatgaacagaatgcaaaaagtggtcgtacttcttgcgaacttcaagaagcagcgaatacctgtcgtcttaaccatttcagccttacagctagaagaattcagcaaaactacttttgcgtcttggggaaactctcttgccttttgggggaacttccattta
cagcagtaaaactcatctttggaagaaagcaacctgtgtactgttggaaaccaagaagtcagcgtcatcaacaacaccctcttttgcgtgaaaatggaaacaaatgtcgttgcgttgggaagaaftctaccacagatataagaataaactctgtgtcatccaagaagataacagattgtatc
ctgctattgtatcatctccagttgggaagtgaatacgtctgcaagttgttaaatatgaagattcaaaactgcagatcttcagcttcaacacgacataaaaactgtgcacttcacttgaagtgtgaagacagattctcatgatcagctgaataacaaaggaaactgaaaacacaagaagaacactt
caagaagctgcataaagccggacatattctatcacgcaagaagtgcaaatgtccctcagctgtgtgggttcagcaaatctgttgcgaagactgttcgcctgaacttcttctgactcgaagattttttctgttaa

More info: Please contact us for any question: contact@proteogenix.fr