

POPDC3-popeye domain containing 3 Gene

Reference: PTXBC022323

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	POP3; bA355M14.1; popeye domain-containing protein 3; popeye protein 3; popeye domain containing 3		

For research use only.

GENE DESCRIPTION	popeye domain containing 3
GENE ID	64208
ORIGIN SPECIES	Human
ACCESSIONS	BC022323
VECTOR	pDONR223
NCBI SYMBOL	POPDC3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaagaattcaagtttatggaagaacctaatagatgaacacccagctgcacaacctggaagcaagagccgaaggagccattatcatctgccagattttatgttagtagtttcatgggtggcagtggaattctcgggctcctttatgtctcagtttgcgggttgggtttctctgtt
ctgctgtctgggctgggttagatgtctgtgacagctgacataatttctggaattttgactgtttgtcatctgcttcacattgttcatattgcatatcaagttcgcagcataacctttgccgagaattccaagtgtgtacagctcctttccagccctgggagctctttgctgtctcagaacg
attgctttgagctctgaagtgttactttgaaaaggaaactgttatgccatgcaggggaaaacttccattgataaactctcctgtctgtttcaggaagatcagagtgacagttgatggcgaattctgcattacattttccctcagttctggaattctcctgagtgagattcactgagacca
cagagggaagcattttcaggaacctcactgcagaaactgattgtcgaatgtgtctgtgaggagaaagaattatctgctcttgcctcagcaltcgtacatctccgccttttccagtgctaaattggcagtgacattgcagataaactctatgccttgaatgacagggtatataaggaaaa
gatatcactatgatattcggctacccaacttctatcaaatgtcaactccagaaatcgcagatcacccctgacacaacattttcagaattccagacgatactgtgataaatga