

**KIAA1598-KIAA1598 Gene**

Reference: PTXBC022348

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms:	KIAA1598; shootin-1; shootin1; shootin 1		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	KIAA1598
<b>GENE ID</b>	57698
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC022348
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	KIAA1598
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

**ORF (Only for reference, product provided without STOP codon )**

atgccagtgttggtctgtatccagtgtaacaaaaacagccttgaacaagaaaactctggagcagaattcaacagcccgtcccccaacactgagccaggtgaagggccccgtaattggaaggatgcacaagttccaagttacgtttcagcctccagtagcattggatgcagga  
aaaaatacattgacggtgaaaaacaagccgaaccagtttagttttatctctgtttctacacatgaaccccaaccaagaccaggtgctgaaaaagatccaactcaacacaagggatgaaggcgaattcaaccgaaaaacaagaagacagcattgaaaacgtgagagagaca  
gacagctccaactgctga