

METTL6-methyltransferase like 6 Gene

Reference: PTXBC022400

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	methyltransferase-like protein 6; methyltransferase like 6		
For research use only.			

GENE DESCRIPTION	methyltransferase like 6
GENE ID	131965
ORIGIN SPECIES	Human
ACCESSIONS	BC022400
VECTOR	pDONR223
NCBI SYMBOL	METTL6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctcttgcagggaaggcgtcaggcaaggattcaccctgaagaaggaggaaactgaaagagaccaacttgggtgtctgattttaaacagcagaaattggaacaaggctcagaaaaattgggatctttttacaaaagaatagcactaattcttcaagacagacact
ggaccaccagagaggttggaggagctaaagatcatgtagagaggttgaagatcaaaagttaacaatgcttgaagctggctgtgggttggaaactgtttatccacttttagaagaagatccgaatcttgcctatgcctgtgattttctccaagagccattgaatatgttaagcaaaatctttata
tgatacagaagatgcaagggtattccagtgtgactgactaaagatgattcttggatcatgtaccgccagagctgtggtggtgtgtatgttgatattgtgtgtcagctgttcacccgataagatgcacctgtcttacaacatttacaagggtattaaaaccaggcaaaagtgtcttgcgtg
actacggactgtatgatcatgccatgcttaggtttaagccagcagcaaaactggagaaaactttatgttagacaagatgggaccagatcataatttttactgatgacttctggtcagctctttatggacacagggtatgaagaagtgttaaacgagtatgtttcagagacggtgaataaa
aaagaaggcctgtga