

SH2D1B-SH2 domain containing 1B Gene

Reference: PTXBC022407

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EAT2; SH2 domain-containing protein 1B; EAT-2; EWS/FLI1-activated transcript 2; SH2 domain-containing molecule EAT2; SH2 domain containing 1B		

For research use only.

GENE DESCRIPTION	SH2 domain containing 1B
GENE ID	117157
ORIGIN SPECIES	Human
ACCESSIONS	BC022407
VECTOR	pDONR223
NCBI SYMBOL	SH2D1B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatctgccttaccatggacgtctgaccaagcaagactgtgagaccttctgctcaaggaaggggtggatggcaactttttaagagacagcagtcgataccaggagtcctgtgctctgtctctgtttaaaaatattgtctacataccgaatcttcagagagaacacgggta
ttacaggatagacagctcagaagggtctccaaaacaggtcttccaagcctaaaggaaactgatctccaaatttgaacccaatcaggggatgggttcaccttttaagccaataaagagaaccagccccagcttgagatggagaggattgaaattagagttggaacatttgtgaacagt
aacagcgattatgtggatgtcttgccttaa