

MOBP-myelin-associated oligodendrocyte basic protein Gene

Reference: PTXBC022471

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	myelin-associated oligodendrocyte basic protein		
For research use only.			

GENE DESCRIPTION	myelin-associated oligodendrocyte basic protein
GENE ID	4336
ORIGIN SPECIES	Human
ACCESSIONS	BC022471
VECTOR	pDONR223
NCBI SYMBOL	MOBP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtcagaaacggccaaggagggtccagactctccaaaaccagaagtactccgaacttcagcatacactgctgcccgccttcaccttctcaattccaagaaggagatagtgatcggaaatacagcatctgaagagcggctgcttaccagaagaagaggaggactg
gatctgctgcctgccagaagaccagattgaaaaggaagatcaggccaacccaagaagaagtga