



Reference: PTXBC022478

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HsT17563; Meflin; immunoglobulin superfamily containing leucine-rich repeat protein; mesenchymal stromal cell- and fibroblast-expressing Linx paralogue; immunoglobulin superfamily containing leucine rich repeat		

For research use only.

GENE DESCRIPTION	immunoglobulin superfamily containing leucine-rich repeat
GENE ID	3671
ORIGIN SPECIES	Human
ACCESSIONS	BC022478
VECTOR	pDONR223
NCBI SYMBOL	ISLR
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]