

GNG4-guanine nucleotide binding protein (G protein), gamma 4 Gene

Reference: PTXBC022485

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	guanine nucleotide-binding protein G(I)/G(S)/G(O) subunit gamma-4; guanine nucleotide binding protein (G protein), gamma 4; guanine nucleotide binding protein 4; G protein subunit gamma 4		

For research use only.

GENE DESCRIPTION	guanine nucleotide binding protein (G protein), gamma 4
GENE ID	2786
ORIGIN SPECIES	Human
ACCESSIONS	BC022485
VECTOR	pDONR223
NCBI SYMBOL	GNG4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaagaggcatgtctataacagcaccactagcatctcccaagccaggaaagctgtggagcagctaaagatggaagcctgtatggacagggtcaaggtctccaggcagctgcggacctctggcctactgtgaagctcagctgcgggaagatcctcatcattccagtgcctgca
tcagaaaacccctttcgcgagaagaagttctttgtaccattctctaa