



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

PRNP-prion protein Gene

Reference: PTXBC022532

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

ASCR; AltPrP; CD230; CJD; GSS; KURU; PRIP; PrP; PrP27-30; PrP33-35C; PrPc; p27-30; major prion protein; alternative prion protein; CD230 antigen; prion-related protein

For research use only.

GENE DESCRIPTION	prion protein
GENE ID	5621
ORIGIN SPECIES	Human
ACCESSIONS	BC022532
VECTOR	pDONR223
NCBI SYMBOL	PRNP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgaacctctgctcgtgatgctggtctcttctgtggccacatggagtgactctgggacctctgcaagaagcgcgaagcctctggaggatggaacactctggggcgacgcgataccggggcaggcgacgcctctggaggcaacctctaccacctcaggcggtgtgtgctggggca
gctcatggtgtggctggggcgagcctcatggtgtgtgtgtggggcgagcccatgtgtgtgctggggacagcctcatgtgtgtgctgggtctgcaaggaggtggcaccacagtctggtgaacaagccgagtgaagccaaaacacatgaagcacatgctgtgtctgcagca
ctggggcgagatggtggggggccttggcgctacgtctgggaagtgccatgacagggcccatcatataattcggcagtgatctaggacaccttactatctgtaaacatgcaccgttaccaccaacagtgtactacagcccatgatgatgacacacagacaactgtgtgac
gactctgcataatcaatcaagcgcacacgctgccaccacaacccaaaggggagaaacttaccgagaccagcttaagatgatggagcgcgtgtgtgagcagctgtatcaccagctacgagaggggaatctacagcctattacaagagagatcgagcatgtctctctctct
ccacctctgactcctcctgactctctcctcatctcctgactcctgagatga

More info: Please contact us for any question: contact@proteogenix.fr