

WBP2NL-WBP2 N-terminal like Gene

Reference: PTXBC022546

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GRAMD7; PAWP; postacrosomal sheath WW domain-binding protein; WW domain-binding protein 2-like; WBP2 N-terminal like		

For research use only.

GENE DESCRIPTION	WBP2 N-terminal like
GENE ID	164684
ORIGIN SPECIES	Human
ACCESSIONS	BC022546
VECTOR	pDONR223
NCBI SYMBOL	WBP2NL
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggatgacagccacacgagaaacggcgcgagccctatccctaacgggtgaaggtctctgaagcgggtctccgaatgtgagctctctccacagcagatcagaaggtcaaatgtcttttagtgtagaaagacaggaacattgtttctcattcataccgggtgattttcataa
cttcattgctccatcagtgatcccatgttgccttttatgatgccatttgatctgacgaacctcactgttgacaacacagtttgcgcaaacctcattaagggaactattcagcagctccatgtgtgctgggaaggacaagctacttttaaatagttctcagaatggagtgccattgaatt
gcccagttgatggfgaagctgcctctgctgctgcccaggatttccacttagaaccttaaatgactgggtcagctctatgggaatttatgtaattactggggaagggaatgtgcactccacagatgcttgcattgtctatggagccccacctgcaggataggagccccacctccc
ggatacggagccccacctgcaggataggagcccaacccgtaggaaatgaaggcccgctgtgggatacagagcctcactgtgcgatalggagccccacctctggatacggagccccacctcaggataggagccccacctctggataggaa
ccccacctctcgataggagccccacctctcgataggagccccacctgcaggaaatgaaggcccgctgcgggatacagagcctcactgtggtcagcaggccagggcctcacgaatctacagcagccaggtctctgaaacgaggctctctctctctctcctctcaggt
ccattcttaa