

MRFAP1-Mof4 family associated protein 1 Gene

Reference: PTXBC022797

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PGR1; MORF4 family-associated protein 1; Mof4 family associated protein 1; T-cell activation protein; protein associated with MRG of 14 kDa; Morf4 family associated protein 1		

For research use only.

GENE DESCRIPTION	Mof4 family associated protein 1
GENE ID	93621
ORIGIN SPECIES	Human
ACCESSIONS	BC022797
VECTOR	pDONR223
NCBI SYMBOL	MRFAP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcggcccttgacatcgtcagctggcgaaccggaggaagtggaggtgctggagcccgaggagatttcgagcagtttctgctcccggtcatcaacgagatgcgcgaggacatcgcgtcgtgacgcgcgagcacggcggtgtacgtgcgaaccggagcaagctgtgg
gagatggacaatatgctcatccagatcaaaacgcaggtggagcctcgaggagagcgcctcaaccacctccagaacccgggcgacgcgcgagggccggcgccaagaggtgcgagaaggccgaggagaaggccaaggagattgcgaagatggcagagatgctgtg
gagctggtccggcgatagagaagagcagtcgtcgtga